



Being sure from lab to production


Your partner for container closure integrity testing and automated visual inspection.

WILCO AG is a Swiss solution provider of leading inspection systems for various packaging types serving the pharmaceutical industry worldwide for more than 50 years.



Enhanced Leak Detection for Rigid Packaging

Durable, robustly built machines enable you to achieve the best inspection quality for your leak tests on rigid packaging by employing sensitive and dependable methods that adhere to global regulations.



Precise Leak Detection of Flexible Packaging

With accurate and adaptable leak detection technologies for flexible packaging, you can ensure product integrity and compliance, while minimizing rejects and increasing throughput.



Automated Visual Inspection

Elevate operational reliability and reduce false reject rates with our automated visual inspection systems tailored for processing small to large production batches.




Multi-Inspection Solutions

Streamline production lines with custom-designed inspection platforms merging visual inspection and leak detection technologies, maximizing product safety while minimizing equipment footprint.



Automated Media Fill Inspection

Efficiently replace manual turbidity checks with our Media Fill solutions to boost accuracy, productivity and regulatory compliance, ensuring adherence to the highest standards.



360° Services

WILCO's extensive services ensure seamless installation, peak performance, and ongoing operational excellence designed to extend your product's lifecycle to a maximum.


Being sure about precision

Whether for laboratory testing, semi-automatic inspection, or high-speed inline testing, our customized strategies are tailored to meet your specific needs. Partner with us to analyze your manufacturing process and create a personalized inspection solution that aligns with your requirements. Drawing from a diverse portfolio of technologies, we combine them to develop the most effective approaches possible.



Leak Detection

- Differential Pressure
- Headspace Analysis
- Force Sensor
- Mass Spectroscopy



Visual Inspection

- Camera based Visual Inspection
- X-ray based Visual Inspection
- Image processing



Process Analytical Technology (PAT)

- Headspace Analysis as IPC
- Near Infrared Spectroscopy



Handling and Coding

- Object Transportation
- Robot Handling
- Tray Loading
- Object Coding

Get your inspection to the next level

Reach out to us to start a discussion about your inspection requirements. We look forward to assisting you.



Tailor-made solutions

We take the time to understand the needs of our customers. Feasibility studies, process analysis, and the method development based on them form the foundation for a successful conception and introduction of an inspection system. Our solutions are individually tailored to your specific needs.

Continuous innovation

As a pioneer in leak testing and visual inspection, we make the lives of people a little safer every day. Our scientific approach is constantly driving us towards new innovations and durable solutions so that we can deliver ever more accurate and reliable measurement results with our inspection solutions.

Broad range of technology

Thanks to our extensive portfolio of technologies and methods, we can find the ideal combination for every type of packaging. Continuous optimization and expansion enable unique solutions and guarantee future-proof solutions with maximum functionality.

Company overview

- Headquarter in Wohlen, Switzerland
- US branch in North Brandford, CT
- Founded 1971
- ~ 165 employees
- More than 5'000 systems in place
- >800 customers in >80 countries
- ISO 9001 certified
- Member of Bausch+Ströbel group

WILCO AG

Rigackerstrasse 11
5610 Wohlen
Switzerland

T: +41 56 618 43 43
info@wilco.com
www.wilco.com

Made in Switzerland

Bausch+Ströbel group